

**HMI** Briefing

Andrew Smith, HSSU Manager







## **July 2013**

Crumlin hospital apologises for contamination scare



O Thursday 25 July 2013 22.21



"Hospital has apologised to 18 families who were wrongly identified as being at the centre of a contamination scare over a medical scope"

Track and Trace?



Our Lady's Children's Hospital, Crumlin has apologised to 18 families who were wrongly identified as being at the centre of a contamination scare over a medical scope.









2013

## July THE IRISH TIMES

Helpline for patients with CJD concerns attracts 1,300 calls

Fri, Jul 19, 2013, 15:39

Between 10 and 20 patients operated on with instruments used on patient with disease

"The HSE will contact up to 20 patients who were operated on using instruments which had been used on a patient who has been diagnosed with CJD disease"

Track and Trace?




The HSE will contact up to 20 patients were operated on using instruments which had been used on a patient who has been diagnosed at Beaumont Hospital with Creutzfeldt-Jakob (CJD) disease Photograph: Dara Mac Dónaill/The Itish Times



















#### **Before**



Manual Check

#### **After**

| j Washer Contents • □ 👂 🗵 |           |                                  |          |              |          |           |  |
|---------------------------|-----------|----------------------------------|----------|--------------|----------|-----------|--|
| Barcode                   | User Code | Description                      | Quantity | Date In      | Time In  | Loaded By |  |
| v Washer 1                |           |                                  |          |              |          |           |  |
| 6000010                   | N/A       | CONTAINERS                       |          | 1 11/04/2012 | 14:00:57 | AGNES     |  |
| • Washer : 3              |           |                                  |          |              |          |           |  |
| 100418500009              | 500       | RHINO & SWR SET (PLASTICS)-02    |          | 1 11/04/2012 | 14:19:54 | PAUL R    |  |
| 100377400043              | 59        | CABG SET-05                      |          | 1 11/04/2012 | 14:19:36 | PAUL R    |  |
| 100293700031              | 3262      | BARRE INSTRUMENTS SET            |          | 1 11/04/2012 | 14:17:01 | PAUL R    |  |
| 100387500031              | 3322      | BOOKWALTER RETRACTOR SET 2-02    |          | 1 11/04/2012 | 14:17:01 | PAUL R    |  |
| 100302000011              | 3525      | FINE STILLE OSTEOTOMES (TH 38.4) |          | 1 11/04/2012 | 14:16:58 | PAUL R    |  |
| 100368000029              | 3099      | GYNAE ONCOLOGY SET-05            |          | 1 11/04/2012 | 14:15:52 | PAUL R    |  |
| 100287700046              | 517       | ABDOMINAL HYSTERECTOMY SET-01    |          | 1 11/04/2012 | 14:15:48 | PAUL R    |  |
| 100322800049              | 3520      | OPEN UROLOGY SET                 |          | 1 11/04/2012 | 14:15:45 | PAUL R    |  |

#### Information available:

- Electronically & Post-event
- Dataloggers linked to tracking system





## Instrument set lists




#### **Before: Paper based**



## After: Scan of barcode and set list automatically printed





#### **After**

- Electronic (list printed when barcode is scanned)
- No need to search in cabinet
- No confusion on document version






## Process log



#### **Before**



Paper based

#### **After**

| To: 05/04    | /2012     |             |            |          |                                      |
|--------------|-----------|-------------|------------|----------|--------------------------------------|
| 01:Packi     | ng Detail |             |            |          |                                      |
| Unknowr      | п Туре    |             |            |          |                                      |
| Barcode      | UserCode  | Operator    | Date In    | Time In  | Description                          |
| 100398700044 | 531       | EILEEN H    | 05/04/2012 | 08:13:31 | GYNAE LAPAROSCOPY SET-01             |
| 100458400011 | 2772      | LINDA R     | 05/04/2012 | 08:54:05 | 30° LENS 4MM WOLF 433424 DSC         |
| 100457900012 | 2130      | LINDA R     | 05/04/2012 | 09:01:13 | 0° LENS 4MM WOLF 427581 DSC          |
| 100458500018 | 2781      | LINDA R     | 05/04/2012 | 09:05:31 | 30° LENS 4MM WOLF 433413 D5C         |
| 100460500024 | 4030      | LINDA R     | 05/04/2012 | 09:08:52 | HYSTEROSCOPE & SHEATH 5000193167 DSC |
| 100303600088 | 2609      | EILEEN H    | 05/04/2012 | 10:41:07 | GENERAL LAPAROSCOPY SET NEW-01       |
| 100455100026 | 3768      | LINDA R     | 05/04/2012 | 10:44:50 | 0° LENS 10MM WOLF 5000171714 DSC     |
| 100375400073 | 2938      | PAUL R      | 05/04/2012 | 11:42:50 | GENERAL LAPAROSCOPY SET NEW-06       |
| 100452500029 | 3195      | ANTHONY C   | 05/04/2012 | 12:11:05 | 0° LENS 5.3MM WOLF 468796            |
| 100454900009 | 3177      | LINDA R     | 05/04/2012 | 14:15:57 | 30° LENS 10MM WOLF 436023 DSC        |
| 100475700020 | 6024      | LINDA R     | 05/04/2012 | 14:16:29 | 0° LENS 5MM STORZ 1127861 DSC        |
| 100375200071 | 2936      | DOROTA S    | 05/04/2012 | 14:31:34 | GENERAL LAPAROSCOPY SET NEW-04       |
| 100321000006 | 2393      | JOSEPHINE D | 05/04/2012 | 17:12:41 | NEPHROSCOPE                          |
| 100453000042 | 3194      | DEIRDRE C   | 05/04/2012 | 17:29:14 | 30° LENS 10MM WOLF 457074            |
| 100452400009 | 3304      | DEIRDRE C   | 05/04/2012 | 17:33:15 | 0° LENS 5MM STORZ 930600             |
| 100416900029 | 3889      | JOSEPHINE D | 05/04/2012 | 17:36:06 | URETEROSCOPE STRAIGHT TH 10 & 11-04  |
| 100453200048 | 3677      | DEIRDRE C   | 05/04/2012 | 17:38:59 | 30° LENS 10MM WOLF 5000145118        |
| 100328200015 | 879       | JOSEPHINE D | 05/04/2012 | 17:42:33 | WOLF LITHOCLAST FULL SET             |
| 100465500034 | 6000      | DEIRDRE C   | 05/04/2012 | 17:44:08 | 30° LENS 4MM ACMI 739453             |
| 100453100047 | 3676      | DEIRDRE C   | 05/04/2012 | 17:48:29 | 30° LENS 10MM WOLF 5000145117        |

Number Of Items: 20

- Electronic
- Post-event can be retrieved at the touch of a button







#### **Document List Report**

| Document Numb   | er: T200242                             | Patient ID: 1087474   | Theatre: THEATRE 2 |
|---|---|---|--------------------|
| Document Created By<br>%JESSIE S                                    | Date & Time Created 02/02/2012 14:14:23 | Amendments  |                    |
| Code: 1003567<br>Lot No: 00035<br>Quantity: 1<br>Operator: %JESSIE: | Item:<br>Supplier:                      | TRACKER<br>BASIC CARDIAC SET-04<br>STERILE SERVICES           |                    |
| Code: 1003774<br>Lot No: 00033<br>Quantity: 1<br>Operator: %JESSIE: | Item:<br>Supplier:                      | TRACKER<br>CABG SET-05<br>STERILE SERVICES                    |                    |
| Code: 1004160<br>Lot No: 00027<br>Quantity: 1<br>Operator: %JESSIE: | Item:<br>Supplier:                      | TRACKER CARDIAC MAMMARY RETRACTOR-03 STERILE SERVICES         |                    |
| Code: 1003741<br>Lot No: 00049<br>Quantity: 1<br>Operator: %JESSIE: | Item:<br>Supplier:                      | TRACKER<br>STERNAL SAW BATTERY POWERED-01<br>STERILE SERVICES |                    |

Number Of Items:

| Document Numbe                     | r: T200243                              | Patient ID: 1087474           | Theatre: THEATRE 2 |
|------------------------------------|---|-------------------------------|--------------------|
| Document Created By<br>%JENNIFER I | Date & Time Created 02/02/2012 15:29:26 | Amendments                    |                    |
| Code: 1002962                      | Тур                                     | DE: TRACKER                   |                    |
| Lot No: 00001                      | Ite                                     | m: CHAUX RETRACTOR CARDIAC-01 |                    |
| Quantity: 1                        | Suppli                                  | er: STERILE SERVICES          |                    |
| Operator: %JENNIFE                 | R I Remar                               | ks:                           |                    |

Number Of Items:

Total Items:

#### **Electronic**

- Post-event Can be retrieved at the touch of a button
- Linked to instrument tracking system





# Loan Sets Major challenge for traceability



#### The Big Challenge!!


- Shared among hospitals (contents always changing)
- Traceability is very challenging
- Can come from a commercial loan set provider or another hospital (often not barcoded)





### Loan set Checklists







- Paper based (hand written)
- No certainty that list matches tray contents









- Electronic, fully legible
- Up to date
- List can be pulled from MS1 Cloud by scanning the tray/set barcode

|     | SE IF PACKING DAMAGED                    |           |             | ER CODE PAGE  |
|-----|--|-----------|-------------|---------------|
| 100 | 2012-04-72                               | Place A 1 | lick In The | Box As Requir |
| 1   | LCP BENDING PIN                          |           | 1           |               |
| 1   | LCP DRTILL SLEEVE                        |           | 1           |               |
| 1   | SCREWDRIVER SHAFT STARDRIVE & HOLDING SL | EEVE      | 1           | No.           |
| 1   | UNIVERSAL DRILL GUIDE                    |           |             |               |
| 0   | BOTTOM TRAY:                             |           |             |               |
| 1   | SCREWDRIVER SHAFT STARDRIVE SHORT        |           | 1           |               |
| 1   | DRTILL BTT 2.0MM                         |           | 1           |               |
|     | DRILL BIT LONG 1.5MM                     |           | 1           |               |









## Loan set – Best Practice Example (PEI)

### PEI Loan sets with GS1 barcodes



|               | MES'S HOSPITAL DUBLIN  Date To: 22/03/2016 |            |  |          |     |            |  |  |
|---------------|--|------------|--|----------|-----|------------|--|--|
| 22-March-2016 |  |            |  |          |     |            |  |  |
| Time          | Barcode                                    | User Code  | Description                                  |          | Qty | Operator   |  |  |
| 13:01:54      | 100563000001                               | ASREXP07IN | ASREXP07IN SET 7, TRAY 4 OF 4                |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100553300003                               | MORUNC02IN | MORUNCO2IN MORELAND UNCEM SET 2, TRAY 1 INST | RUMENTS  | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564200000                               | SROM01IN   | SET 1, TRAY 5 OF 7, IMPACTOR/EXTRACTOR INSTS |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564900000                               | MOREXT01IN | MORELAND EXTRACTION SET 1, TRAY 1 INSTRUMENT | S.       | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564800000                               | SOLUTN01IN | SOLUTION THINSHAFT REAMERS, TRAY 6 OF 6      |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564700000                               | SOLUTN01IN | SOLUTION INSTRUMENTS TRAY 5 OF 6             |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564500000                               | SOLUTN01IN | SOLUTION TRIAL RIGHT. TRAY 3 OF 6            |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100562800001                               | ASREXP07IN | ASREXP07IN SET 7. TRAY 3 OF 4                |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564300000                               | SOLUTN01IN | SOLUTION FEMORAL REAMERS TRAY 1 OF 6         |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564600000                               | SOLUTN01IN | SOLUTION TRIAL STANDARD. TRAY 4 OF 6         |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100562700001                               | ASREXP07IN | ASREXP07IN SET 7, TRAY 1 OF 4                |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100562900001                               | ASREXP07IN | ASREXP07IN SET 7, TRAY 2 OF 4                |          | 1   | CAROLINE F |  |  |
| 13:01:54      | 100564400000                               | SOLUTN01IN | SOLUTION TRIAL LEFT. TRAY 2 OF 6             |          | 1   | CAROLINE F |  |  |
|               |  |            |  | SubTotal | 13  |            |  |  |
|               |  |            |  | Total    | 13  |            |  |  |


#### Loan sets from PEI scanned into Wash Room

- No need to re-key information (huge time saving)
- Full audit trail





# Barcode (keydot) being added to supplementary instruments










# Overview of the number of procedures year on year in SJH HSSU \*\*Retween 2014 to 2018 increase

\* Between 2014 to 2018 increase in workload with less staff



Anecdotal evidence:
Moved from 5,000
Transcriptions per day
in one hospital





## **Key Benefits**

### National traceability solution for Instrument trays and endoscopes

- Automatic Tracking of Instrument Sets
- Much easier to share Loan Sets
- Tray Checklists printed when scanned
- All records stored digitally (huge paper saving)
- Link between tracking system and theatre to patient (EHR)
- Improved Workflow (time savings)



Reporting and management of set inventory






- Adding Keydot barcode to all supplementary instruments
- Adding additional areas such as the Cath Lab & OPD which use sterile surgical sets
- Requesting all loan set providers to be compliant











#### **Andrew Smith**

HSSU, St James's Hospital James's Street, Dublin, Ireland

Tel + 353 1 416 2557

eMail ASmith@stjames.ie



